





# Autograph Plus



Minimum Fitting Height 14mm, 16mm, 18mm



# Autograph Plus



Fitting Height	A	B	C (Invisible Marking)
14mm	2	7	+14
16mm	2	9	+16
18mm	2.5	11	+18

Fitting Point

Distance Reference Point

Prism Reference Point

Near Reference Point

Ink Marking: Blue  
Invisible Marking: Red

Minimum Fitting Height 14mm, 16mm, 18mm

Material	Invisible Marking
Hard Resin 1.5	15 +16 250
Transitions® 1.5	TR +16 250
Polarized 1.5	PZ +16 250
Drivewear® Activated by Transitions® 1.5	DW +16 250
Trivex™ Lens Material 1.53	TX +16 250
Transitions® 1.53	TTX +16 250
Polarized 1.53	PTX +16 250
Polycarbonate 1.59	PC +16 250
Polycarbonate 1.59 Transitions®	PT +16 250
Polycarbonate 1.59 Polarized	PPZ +16 250
Drivewear® Activated by Transitions® 1.59	PD +16 250
SuperLite™ 1.6	16 +16 250
SuperLite 1.6 Transitions®	T6 +16 250
SuperLite 1.6 Polarized	PZ6 +16 250
SuperLite 1.67	67 +16 250
SuperLite 1.67 Transitions®	T67 +16 250
SuperLite 1.67 Polarized	PZ67 +16 250
SuperLite 1.74	74 +16 250