







Shamir Autograph III[®]



Minimum Fitting Height: 11mm



Shamir Autograph III®



Fitting Point

Distance Reference Point

Prism Reference Point

Near Reference Point

Ink Marking: Blue
Invisible Marking: Red

Minimum Fitting Height: 11mm

Material	Invisible Marking xx	Material	Invisible Marking xx
Hard Resin 1.5	15	Polycarbonate 1.59 Polarized	PPZ
Transitions® 1.5	TR	Drivewear® 1.59	PD
XTRActive™ 1.5	X5	SuperLite™ 1.6	16
Polarized 1.5	PZ	SuperLite 1.6 Transitions®	T16
Drivewear® 1.5	DW	SuperLite 1.6 XTRActive™	X6
DLC™ 1.53	TX	SuperLite 1.6 Polarized	PZ6
DLC 1.53 Transitions®	TTX	SuperLite 1.67	67
DLC 1.53 Polarized	PTX	SuperLite 1.67 Transitions®	T67
Polycarbonate 1.59	PC	SuperLite 1.67 XTRActive™	X67
Polycarbonate 1.59 Transitions®	PT	SuperLite 1.67 Polarized	PZ67
Polycarbonate 1.59 XTRActive™	PX	SuperLite 1.74	74