




## Shamir SV III As Worn™



Fitting Point  
Prism Reference Point



Ink Marking: Blue  
Invisible Marking: Red

Material	D Invisible Marking
1.5 Hard Resin Clear	15
1.5 Polarized	PZ
1.5 Transitions® Signature®	TR
1.5 Transitions® Vantage®	VA
1.5 Transitions® XTRActive®	X5
1.5 Transitions® Drivewear®	DW
1.53 DLC™ Clear	TX
1.53 DLC™ Polarized	PTX
1.53 DLC™ Transitions® Signature®	TTX
1.53 DLC™ Transitions® Vantage®	VX
1.53 DLC™ Transitions® XTRActive®	XX
1.56 Mid Index	56
1.59 Polycarbonate Clear	PC
1.59 Polycarbonate Polarized	PPZ
1.59 Polycarbonate Transitions® Signature®	PT
1.59 Polycarbonate Transitions® Vantage®	PV

Material	D Invisible Marking
1.59 Polycarbonate Transitions® XTRActive®	PX
1.59 Polycarbonate Transitions® Drivewear®	PD
1.6 SuperLite® Clear	16
1.6 SuperLite® Polarized	PZ6
1.6 SuperLite® Transitions® Signature®	T6
1.6 SuperLite® Transitions® Vantage®	V6
1.6 SuperLite® Transitions® XTRActive®	X6
1.67 SuperLite® Clear	67
1.67 SuperLite® Polarized	PZ67
1.67 SuperLite® Transitions® Signature®	T67
1.67 SuperLite® Transitions® Vantage®	V67
1.67 SuperLite® Transitions® XTRActive®	X67
1.74 SuperLite® Clear	74
1.74 SuperLite® Transitions® Signature®	T74
1.74 SuperLite® Transitions® XTRActive®	X74





# Shamir SV III As Worn™





# Shamir Autograph InTouch™



Minimum Fitting Height 15mm, 18mm

